



CloudPAT

Itamar Medical OM2200005

Web Software Operation Manual

WatchPAT® and PAT® are trademarks of Itamar Medical®, Ltd.


When using the marks in our publications and this website, we used the appropriate TM, SM, or [®]. The [®] symbol is used to indicate registrations in the U.S. and the territories listed in <https://www.itamar-medical.com/terms-and-conditions/> and this should not mean that such trademarks are necessarily registered in other countries. For example, WatchPAT® is a trademark of Itamar Medical Ltd., registered in the U.S. and the applicable territories listed in <https://www.itamar-medical.com/terms-and-conditions/> only.



Caution: Federal law restricts this device be sold by the order of a licensed healthcare practitioner only.

Copyright © 2024 by Itamar Medical Ltd.

Software Version: 3.0

This manual and the information contained herein are confidential and are the sole property of  Itamar Medical Ltd. Only Itamar Medical Ltd. or its licensees have the right to use this information. Any unauthorized use, disclosure or reproduction is a direct violation of Itamar Medical's proprietary rights.

DISCLAIMER

Itamar Medical Ltd. shall not be held responsible in any manner for any bodily injury and/or property damage arising from operation or use of this product other than that which adheres strictly to the instructions and safety precautions contained herein and in all supplements hereto and according to the terms of the warranty provided in the License Agreement in Appendix F. Screen layouts are subject to change from time to time without notification and without being altered in this Web software operation manual as long as the intent of the content has not been altered significantly.

This product and/or method of use is covered by one or more of the following US patents: 6319205, 6322515, 6461305, 6488633, 6916289, 6939304, 7374540, as well as any pending US patent applications and corresponding patents and/or applications filed in other countries.



Itamar Medical Ltd.
9 Halamish St., P.O.Box 3579
Caesarea Ind. Park, 3088900, Israel
Tel + 972 4 617 7000
Fax + 972 4 627 5598
www.itamar-medical.com
support@itamar-medical.com



Arazy Group GmbH
The Squire 12, Am Flughafen, 60549 Frankfurt
am Main, Germany

EN ISO 13485:2016

See appendix B for contact information of the regulatory authorized representative

Document Information:

Project Name	CloudPAT Web Software Operation Manual
Document Number	OM2200005
Revision No	9
Applicability Date	July 2025

Revision History:

Date	Revision Level	Author	ECO No.	Change Description
Jan 2021	1			1st Release
Apr 2021	2			GA Release
Sep 2021	3			Version 1.4
Aug 2022	4			Versions 1.5 – 1.7
Mar 2023	5			Versions 1.8 – 3.10
Jan 2024	6			Versions 3.11 – 3.12.5
Jan 2025	7			Versions 3.13 – 3.14.1
May 2025	8			Versions 3.15 – 3.16
July 2025	9			Versions 3.16 – 3.17

Table of Contents

Table of Contents.....	4
1 Introduction to CloudPAT	7
1.1 Intended Use of the CloudPAT (CP) ®	7
1.2 zzzPAT & WatchPAT Interface S/W – Definition.....	7
1.3 CloudPAT S/W – Definition	7
1.4 Overview	7
1.5 Terminology	7
1.6 Minimum Requirements.....	8
2 Setting CloudPAT Access.....	9
2.1 Obtaining an Account	9
2.2 Accessing CloudPAT	9
2.3 Setting [Initial] Password	10
2.3.1 Forgot Password / User Name.....	11
2.3.2 Contact Us.....	12
2.3.3 Email Preferences	14
3 Using CloudPAT as a Provider	15
3.1 Main Access Screen	15
3.2 Dashboard.....	15
3.2.1 Enhanced Patient Care Pathway.....	16
3.2.2 Screening Patient Care Pathway.....	17
3.2.3 Sleep Apnea Severity Scale	18
3.2.4 Central Sleep Apnea bar chart.....	19
3.2.5 Screening Patients STOP BANG bar chart.....	20
3.2.6 Screening Patients ESS bar chart	21
3.2.7 Longest AFib duration pie chart	22
3.3 Managing Patients	22
3.4 Creating a New Patient	25
3.4.1 Creating a New WatchPAT Patient	26
3.4.2 Creating a New Screening Patient	27
3.5 Importing Patients	27
3.6 First Time Installation	28
3.7 Initialize WP200U/WP300 Device.....	31

3.8	Initialize WatchPAT™ ONE Device	33
3.9	Initialize WatchPAT™ ONE Multi-Night Device.....	34
3.10	Initializing a WP400 device	35
3.11	Device Initialization Details.....	36
3.12	Uploading a Study	36
3.13	Normal Upload Mode (WP200U/WP300)	37
3.14	Manual Mode (WP200U/WP300).....	39
3.15	Manual Upload of WatchPAT ONE Study	40
3.16	Searching Studies and Reports	40
3.17	Study Management	41
3.17.1	Understanding the Study Screen	42
3.17.2	Attaching Files to a Patient.....	43
3.17.3	Changing Interpreting Physician.....	44
3.17.4	Interpretation Pools.....	45
3.17.5	View Analysis File.....	46
3.17.6	Request for a Script	47
3.17.7	Changes of Study Parameters.....	48
3.18	Using CloudPAT as an Office Admin	48
	Creating office Providers:	48
	Creating Office Interpreting Physicians:	48
	Creating Office Sleep Physicians:.....	49
3.19	SleePATh Questionnaires.....	49
3.19.1	Learning center:.....	51
3.20	IDTF Customers - HSAT Order Form	51
3.21	Viewing Reports	52
3.22	Managing Devices.....	53
3.22.1	Devices Table	54
3.22.2	Device Details	55
3.23	Managing Referring Physician	55
3.23.1	Referring Physicians Table	56
3.23.2	Referring Physician Details	57
3.24	Managing Lists	57
3.25	SleePATh Communications Tab	58
3.26	Using CloudPAT as a Therapy Provider.....	59

4	Using CloudPAT as an Interpreting Physician	60
4.1	Main Access Screen	60
4.1.1	Studies Assigned to the Interpreter.....	61
4.1.2	Reviewing a Study.....	62
4.1.3	Viewing a Completed Study.....	66
4.2	Manage Templates & NPIs.....	66
4.2.1	Templates	67
4.2.2	NPIs	68
4.2.3	Signature	68
4.2.4	States Eligibility.....	68
4.3	Using CloudPAT as an Interpreting Physician Assistant	69
5	Appendices	70

1 Introduction to CloudPAT

1.1 Intended Use of the CloudPAT (CP) ®

The CloudPAT (CP) is a web-based, cloud system for patient care pathway management specifically for patients suspected to have sleep-related breathing disorders. The CloudPAT is intended to display and store the medical device data from WatchPAT by interfacing with zzzPAT modules. The system allows physicians to enter patient data.

1.2 zzzPAT & WatchPAT Interface S/W – Definition

zzzPAT and WatchPAT Interface are software packages used with the WatchPAT™ device to aid in the diagnosis of sleep related breathing disorders, detecting REM, Light Sleep, Deep Sleep and Wake stages and measuring snoring intensity and body position states. The S/W displays the signals recorded by the WatchPAT™ device, automatically identifies breathing disorder events, records sleep stages; tracks snoring and body position data; and then generates a comprehensive report for the physician.

1.3 CloudPAT S/W – Definition

1. The CloudPAT (CP) is a HIPAA-compliant, web-based cloud application providing a management system for sleep/breathing disordered patients from diagnostic to treatment compliance.
2. CP uses data that has been taken from the WP device. CloudPAT uses the WatchPAT Interface package for the WatchPAT device.
3. CP is intended to display and store the medical device data without controlling or altering its functions or parameters. CloudPAT does not modify, interpret, or add value to the data or the display of the data.
4. CP allows management of the study by different practitioners, including the preparation and uploading of the study by one group of sleep technicians/physicians and the review and interpretation of the sleep study by other sleep specialists.
5. CP is compatible with the WatchPAT mobile application - assigning digital questionnaires, sending notifications, and presenting data collected via SleepPath.
6. CP supports external APIs such as FHIR and Redox interconnectivity.

1.4 Overview

CloudPAT uses the zzzPAT's analysis. In cases where the night's data needs to be viewed and automatically detected events need to be revised manually, the study data can be downloaded to a local PC using the zzzPAT.

This manual provides the information necessary for routine use of the CloudPAT software.

1.5 Terminology

CloudPAT supports various types of users and entities.

- Office – a practice responsible for administering the tests. Office can include one or more referring physicians, as well as technicians and administrative staff.
- Referring Physician – the primary physician or dentist prescribing the sleep test.
- User – typically a Provider such as a technician or a nurse, but might also be a physician. This is the most common type of user using CloudPAT functions such as Creating a new study, Uploading a study, Viewing reports and Managing Patient Information.
- Interpreting Physician– a physician (typically Board-certified) performing the sleep study analysis, interpretation, diagnosis and recommendations. CloudPat functions such as reviewing,

interpreting studies, managing templates for diagnosis and recommendation may also be performed.

- Sleep Physician – a combination of Provider and Interpreter.
- Patient – the person undergoing the sleep study.

1.6 Minimum Requirements

An internet connection using a PC workstation with Microsoft Win10 and above, running at a minimum upload speed of 1 Mbps and supported browser (Edge, Chrome or Firefox) is required.



Note

Figures included in this Operation Manual are for reference only and the actual design and appearance might be different in the final product.

2 Setting CloudPAT Access

2.1 Obtaining an Account

Usage of CloudPAT is intended for WatchPAT customers only (using all supported devices). A secure and individual login account must be assigned by an Itamar administrator.

To obtain a User Name and Password, contact your Sales Representative or send an email to CloudSupport@itamar-medical.com. After a vetting process, you will get an email invitation containing your login information.

CloudPAT supports SSO connectivity. If enabled for your organization, the login shall be done through your SSO portal. Please contact your system administrator for further details.

CloudPAT supports the following languages based on user settings (including emails, zzzPAT viewer, Sleep Report and Exports to CSV): English, French.

2.2 Accessing CloudPAT

All subsequent access to CloudPAT will be according to the server location:

Server location	Direct link
United states	https://cloudpat.itamar-medical.com/
European Union	https://cloudpat-eu.itamar-medical.com/
United Kingdom	https://cloudpat-uk.itamar-medical.com/
Australia	https://cloudpat-au.itamar-medical.com/

Access to CloudPAT begins with the Sign In screen and includes a link to Itamar's Terms and Conditions Screen where User Name and Password are entered.

Sign In Screen

After several failures to sign into CloudPAT, your account may be temporarily locked. If such an event occurs, CloudPAT will send a notification message to your email account.

An office you belong to may restrict CloudPAT access to a specific IP address range, or to a particular country group within their region.

If you do not meet the access settings, your login will be blocked, and the proper error message will be displayed.

To set up login restrictions, please contact CloudPAT Support.

After a successful login, a “What’s New” message may be displayed, listing the new features of the latest version and a banner may be displayed at the top of the screen.

2.3 Setting [Initial] Password

Upon initial account setup, you will be prompted to enter a temporary password.

CloudPAT will prompt you to pick a new password upon first-time sign in.

The minimum password length is eight characters, and the password must contain characters from three of the following four categories: English uppercase characters (A through Z), English lowercase characters (a through z), digits (0 through 9), non-alphabetic characters (for example: ! , \$, # , %).

The new password shall be different from any of the user’s previous passwords.

CloudPAT supports additional security measures, and allows users to log in with MFA through a text message sent to the patient’s mobile phone. This option needs to be enabled by the system’s admin upon request.


Change Password Screen

**Note**

All asterisk (*) tagged fields are mandatory and should be filled.

2.3.1 Forgot Password / User Name


If you happen to forget your password, simply click on the “Forgot your password?” link on the Sign In screen, enter your User Name and an email will be sent to you with a link to reset password.



Forgot Password

Please fill in your username to continue with password recovery. If you don't know your username please contact us.


* USERNAME

 Need help signing in?

Cancel Continue

Forgot Password Screen


If you happen to forget your User Name, simply click on the “Forgot your username?” link on the Sign In screen, enter your Email and an email will be sent to you with the associated User Name.



Forgot Username

Please fill in your email address and associated usernames will be sent to you

* EMAIL

 Need help signing in?

Cancel Continue

Forgot User Name Screen

2.3.2 Contact Us

The Contact Us screen includes two communication channels as well as links for downloading the latest CloudPATio client and WatchPAT 300 Device Drivers.

The screenshot shows the 'Contact Us' dialog box. At the top, there are links: 'The latest CloudPATio client can be downloaded [here](#)' and 'Download the latest WatchPAT 300 device drivers [here](#)'. Below these links is a list of four contact options, each with a radio button and an icon:

- General Inquiry**: Represented by a clipboard icon. This option is selected, indicated by a blue dot in the radio button.
- Send Study for Clinical Support**: Represented by a clipboard icon.
- By Phone**: Represented by a telephone handset icon.
- Request for record removal from CloudPAT server**: Represented by a clipboard icon.

 To the right of the list is a form with the following fields:

- NAME**: A text input field.
- EMAIL**: A text input field.
- OFFICE**: A dropdown menu.
- MESSAGE**: A large text area for the inquiry.

 At the bottom of the dialog are 'Cancel' and 'Send' buttons.

Contact Us Screen – General Inquiry

For administrative support and help with any problems you experience with CloudPAT, please contact us by clicking the “Contact Us” link from any window.

The “General Inquiry” creates a ticket for Itamar-Medical personnel to address.

Please contact customer support by selecting the “By Phone” option if any difficulties arise. This screen contains the up-to-date phone numbers.

This screenshot shows the 'Contact Us' dialog box with the 'By Phone' option selected. The layout is similar to the previous screenshot, but the radio button for 'By Phone' (represented by a telephone handset icon) is now selected with a blue dot. The 'General Inquiry' option is no longer selected. The form fields on the right (NAME, EMAIL, OFFICE, MESSAGE) are still present but empty. The phone number 'US and Canada: TEL: 1 888 748 2627' is displayed above the 'By Phone' option. At the bottom, the buttons are 'Cancel' and 'OK'.

Contact Us Screen – By Phone

If enabled, the “Contact Us” menu has the option to send a study to consultation with ZOLL’s clinical support. This option is available only if the user is viewing the “Patient Details” screen, and the patient has at least one uploaded study when selecting the Contact Us menu option.

If enabled, the “Contact Us” menu has the option to request for record removal from CloudPAT servers. Once this office setting is turned on, this option will allow users to specify the desired criteria for the removal.

2.3.3 Email Preferences

There is a link called “Manage User Preferences” under the User menu which opens a designated screen listing the email options for the user. The user may customize the options. Email options are role based and not identical to all system users.

User Email Preferences

X

toggle off all options that you prefer not receive email notifications for

There is a new study awaiting interpretation	<input checked="" type="checkbox"/>	Indication about a study that has been assigned for your interpretation.
The report has been unlocked	<input checked="" type="checkbox"/>	Indication about a study that has been locked by you has been unlocked.
The study has been reassigned	<input checked="" type="checkbox"/>	Indication about a study that has been assigned for your interpretation is now re-assigned.
There is a new study awaiting pool interpretation	<input checked="" type="checkbox"/>	Indication about a study that has been assigned for a pool.
Daily Interpretation Pool Summary	<input checked="" type="checkbox"/>	Periodic summary email in relations to pending pool studies.
A new study has been uploaded	<input checked="" type="checkbox"/>	Indication that a study you have initialized has been uploaded to CloudPAT.
Analysis has failed	<input checked="" type="checkbox"/>	Indication that a study you have initialized has failed to run analysis.
A new report has been completed	<input checked="" type="checkbox"/>	Indication that a study you have initialized has been interpreted and locked.
Screening Uploaded	<input type="checkbox"/>	Indication that a screening study you have initialized has been uploaded to CloudPAT.

Cancel

Save

Email Preferences

3 Using CloudPAT as a Provider

3.1 Main Access Screen

Here are the main actions that can be selected from this screen:

- Dashboard - monitors patient progress and gives physicians and staff better control over patient care.
- New Patient – prepares a new patient & initializes WatchPAT™ device for a new sleep study.
- Upload Study – uploads sleep study files from WatchPAT™ device to the CloudPAT web server and subsequent redirection (e.g., interpretation by an Interpreting Physician) according to pre-determined instructions.
- Patients – enables the user to see which sleep studies have been performed on patients, alter the patient's information and/or generate a follow-up sleep study for an existing patient.
- Devices – Displays a table of the associated devices. Also allows access to more details for a particular device.
- Managing Referring Physicians – allows the user to directly manage the referring physicians associated with the office.
- Configuring Managed Lists - allows the user to manage the following lists: Medications, Diseases, Statuses, and Insurers. The user may also search for existing entries from global and local office level lists or add new ones.

3.2 Dashboard

The Dashboard includes the following items:

3.2.1 Enhanced Patient Care Pathway

The Enhanced Patient Care Pathway is a stacked bar chart displaying statistics from the selected patient group. It follows the patient through four stages, beginning at check-in until assessment of their Sleep Apnea.

At each stage along the pathway, the total number of patients is listed, along with the percentage from the previous stage.

Each stage is comprised of 1-3 sections that, when combined, equals to the number of patients in the first section of the previous stage.

Selecting any of the sections will display a list of patients based on the selected criteria, and selecting a patient will display the Patient Details Screen.

Tooltip is displayed next to the segment when the mouse hovers over the icon.

The results can be filtered by Office, Referring Physician and Last Status Update.


The following stages are displayed in the chart:

1. Patients checked in: The number of patients in selected group (percentage is always 100%).
2. WP assigned:
 - a) WP assigned (main): Patients assigned with initialized WatchPAT device, including patients who have returned the device, and the WP study was uploaded to CloudPAT but not yet analyzed.
 - b) WP not assigned: Patients not yet assigned a WatchPAT device.
3. Tested with WP:
 - a) Tested with WP (main): Patients that have successfully completed the WatchPAT Study and the data was successfully analyzed.
 - b) WP Study Failed: Patients' WP study was uploaded to CloudPAT, but no AHI was calculated.
 - c) Not tested yet: Patients' WP study not yet uploaded to CloudPAT.
4. Sleep Apnea:
 - a) Apnea Positive (main): Patients positively diagnosed with sleep apnea, i.e., AHI above the office-defined threshold.
 - b) Apnea Negative: Patients with AHI below the threshold.
 - c) Report Not Ready: Patients with successfully uploaded studies to CloudPAT but the report is not ready.

Patients Care Pathway

AHI Threshold used for calculation: 12

● Patients Created ● WP Assigned ● WP Not Assigned ● Tested With WP ● Failed WP Study ● Not Tested with WP ● Sleep Apnea ● No Sleep Apnea ● Report Not Ready




Enhanced Patient Care Pathway

3.2.2 Screening Patient Care Pathway

If your office is configured for screening questionnaires, you can create a new screening patient by compiling the answers to a series of questions from the WatchPAT mobile app.

CP provides a screening funnel widget based on stacked rows displaying the patients' progression in the screening process. Each bar represents the status in respect to the relevant stage.

1. Screening Created
2. Screening Assigned
3. Screening Loaded
4. WatchPAT Decision
 - a) Assign WP
 - b) No Sleep Study
 - c) Pending Decision.



Screening Patient Care Pathway

3.2.3 Sleep Apnea Severity Scale


The Sleep Apnea Severity Scale is a pie chart indicating the Sleep Apnea Severity Scale for Tested with WP patients containing four levels (according to the office settings):

1. No Apnea
2. Mild Apnea
3. Moderate Apnea
4. Severe Apnea

When a section is selected, the respective query in the patients' table and the Sleep Apnea filter will display the four levels.

Sleep Apnea Severity Scale

AHI Threshold used for calculation: 5, 15, 30




Sleep Apnea Severity Scale

3.2.4 Central Sleep Apnea bar chart

The Central Sleep Apnea bar chart shown to the right of the Sleep Apnea Severity Scale including these ranges of % AHlc / AHI. The calculations are based on the Central office settings. The thresholds are displayed if they are identical for the requested offices.

Central Sleep Apnea

The specific SA Threshold of each of the participating offices will be used for the calculation AHlc as a percentage of total AHI




Central Sleep Apnea bar chart

3.2.5 Screening Patients STOP BANG bar chart

The Screening Patients STOP BANG bar chart shows information based on data received from SleepPATH screening, within these ranges: 0 – 2, 3 – 4, 5 – 8.

Screening Patients STOP BANG

*Includes only the patients in screening stage without a WP device initialization process.




Screening Patients STOP BANG chart

3.2.6 Screening Patients ESS bar chart

The Screening Patients ESS bar chart shows information based on data received from SleepPATH screening within these ranges: 0 - 5, 6 - 10, 11 - 12, 13 - 15, 16 – 24.

Screening Patients ESS

*Includes only the patients in screening stage without a WP device initialization process.



Screening Patients ESS bar chart

3.2.7 Longest AFib duration pie chart

Disclaimer: The WatchPAT is not intended to be used as a diagnostic device for any cardiac arrhythmia and is not intended to replace traditional methods of diagnosis for cardiac arrhythmias.

The Longest AFib duration is a pie chart indicating AFib episodes duration with WP patients according to four levels:


1. Not Enabled
2. Not Detected
3. Less than one min
4. Greater than one min

When a section is selected, the respective query in the patients table will display.

In cases where the longest AFib episode detected is shorter than 60 seconds, review of the peripheral arterial tonometry signal for an irregular rhythm in the location of the episode is recommended. The WatchPAT may not detect short AFib episodes (<60 seconds).

Longest AFib Duration in Sleep

*The WatchPAT® is not intended to be used as a diagnostic device for any cardiac arrhythmia and is not intended to replace traditional methods of diagnosis for cardiac arrhythmias.



Longest AFib duration pie chart

3.3 Managing Patients

You may browse or edit existing patient details by pressing “Patients” from the main screen which will display the Search for Patients Screen. Click “Search” to see a list of all patients.

You can enter filters to limit the search results.

The search fields are dynamic and so are the table columns.

You may click on the Export button for various csv export options of the list of patients to your PC.

There is an option to toggle between single (latest) study per patient to multiple (show all studies).

You may choose "Select all" and "Un-Select all" in some of the dropdown search fields.

PATIENT ID	NAME	DATE OF BIRTH	OFFICE	PATIENT CREATED DATE	STUDY DATE	STATUS
ronen111			Ronen Office	08-14-2022	08-14-2022	
ronen special ch	000j3kac4444880ac		Ronen Office	08-11-2022	08-11-2022	
ronen_wg300			Ronen Office	08-04-2022		
ronen102			Ronen Office	04-06-2022		
ronen101			Ronen Office	04-05-2022	04-05-2022	
ronen100			office	04-04-2022	04-04-2022	
ronen12	LN, FN		office	04-30-2021		
ronen005			office	04-27-2021		
ronen003	123456789012345, 123456789012345		office	04-20-2021		
ronen002	Watson, Emma	04-06-2009	office	04-11-2021		

Patient Table Screen

Clicking on a row opens the patient's details screen.

If the office enables Screening Patients, there is a Screening Patients Table, which does not include any of the study related search fields and columns.

a Screening Status as search field and column can be added.

There are additional search fields called STOP BANG and ESS, with the same values from the dashboard.

PERSONAL DETAILS	CONTACT INFORMATION	INSURANCE INFORMATION
HEIGHT(IN) 71	STREET street 1	INSURANCE PROVIDER -
WEIGHT(LB) 177	CITY city 1	GROUP NUMBER -
BMI 24.69	ZIP CODE zip 1	OTHER INSURANCE -
NECK CIRC(IN) 16	STATE Alaska	INSURANCE ID NUMBER -
EPWORTH SCORE 18	COUNTRY -	PRE-AUTH NEEDED No
STOP BANG SCORE 5	HOME PHONE -	PRE-AUTH NUMBER -
PACEMAKER No	ADDITIONAL PHONE NUMBER -	PATIENT'S FINANCIAL RESPONSIBILITY -
PREFERRED PATIENT LANGUAGE	NOTIFICATION No	BILLING COMMENT -

Patient Details Screen

Subject to office setting, a Provider user has the ability to delete patients from the CP office using a designated “Remove Patient” button. A confirmation message will be displayed before the patient is deleted.

If the patient has an initialized WP1 device, a confirmation message will be displayed, and once confirmed an automatic unregister operation will occur.

The screenshot shows a patient profile page with a progress bar at the top indicating five stages: 1. Created (10-09-2022), 2. In Testing (initiated) (10-09-2022), 3. Study Loaded (10-09-2022), 4. Sent to Interpretation (03-21-2024), and 5. Study Completed. Below the progress bar, patient details are listed: PATIENT ID, NAME, OFFICE, GENDER (Male), DATE OF BIRTH, AGE, REFERRING PHYSICIAN, COUNTRY CODE, MOBILE PHONE, and EMAIL. At the bottom right, there are two buttons: 'Initialize Device' (green) and 'Remove Patient' (red). Below the patient details, there is a tabbed interface with options: 'Studies & Initializations', 'Details', 'WPT', 'Diseases & Medications', 'Attached Files', 'SleepPATH™', and 'HSAT Medical Order Form'.

Remove Patient button

Providers can change a Patient's ID for uploaded studies only once, with the following conditions:

- The patient has one study
- There are no manual events
- The study is not in in-testing stage

Providers can disable all SleepPATH services from the Patient Card, in the SleepPATH tab.

In case of manual events, if the analysis has been edited or changed, an icon will appear on the patients table with the description of: "The automatic analysis events or stages have been edited."

The screenshot shows the CloudPAT Patients table. The table has columns: PATIENT ID, NAME, DATE OF BIRTH, INSURANCE PROVIDER, OFFICE, REFERRING PHYSICIAN, INTERPRETING PHYSICIAN, PATIENT CREATED DATE, PATIENT SYSTEM ID, and STATUS. A red square icon is visible in the STATUS column for a patient with a date of birth of 07-08-1999. The table is filtered by 'Last Status Update' and 'Two Weeks'. The left sidebar contains search filters for Patient ID, First Name, Last Name, Patient System ID, Post Education Eligibility, Office, WPT Partner, Status, and Interpreting Physician.

Manual Events Indicator Icon

3.4 Creating a New Patient

When creating a new patient, erroneous fields and empty mandatory fields will be highlighted in red.

It is possible to configure specific fields as additional mandatory fields by contacting the system administrator.

The user may attach files which will be associated to the patient's record.

You can save the patient record at any stage and return to it later from the "Patients" screen.

Filling in the patient's name is recommended for easier patient management.

There is an office setting to configure key patient details fields as mandatory or not.

When a new patient is created or edited for an office, saving would not be possible if fields defined as mandatory are left blank.

Patient Details screen includes the following tabs: Details, Diseases & Medications, Attached Files and options for Studies & Initializations, App Questionnaires and Compliance Data

Patient Details tab includes the following sections: Header Section, Personal Details, Logistic Comments, Contact Information, Study Details, Insurance Information, Additional Information, Status Information.

Changes to Patient Details are captured in the Audit Trail Screen by clicking on the Log button.

3.4.1 Creating a New WatchPAT Patient

In order to prepare the WatchPAT™ device for a new patient please select “New WP Patient”.

This will create a new patient record and open the patient screen.

The following steps are involved in preparing for a new study:

- The mandatory fields are Patient ID and the selection of an office, if the user is associated with more than one office. Other mandatory field may be defined at the office level.
- The Patient ID cannot be changed once a device was initiated.
- Pacemaker – select this option if the patient has an implantable pacemaker. CloudPAT will automatically detect segments where the PAT pulses seem to be paced and exclude these segments from the analysis.



Note

Some types of Pacemakers are excluded from the use of the WatchPAT device. See the device's Operating Manual for more details.

Once the patient definition is completed, press “Save and Initialize Device”.

WP Patient Card

3.4.2 Creating a New Screening Patient

- Similarly to WP Patients, if your office is configured for screening questionnaires, it is possible to create a new screening patient by using the New Screening Patient button.
- This would allow getting the patient's answers to a series of questions through the WatchPAT mobile app to be added as a PDF report in the App Questionnaires tab.
- A screening patient has additional mandatory fields to be used to identify in the WatchPAT mobile app such as Email/DOB (per office settings) and mobile phone number.
- Most of the patient tabs are similar to a WP patient card except Compliance Data and Study & Initializations tabs.
- Clicking on the Initialize Screening button triggers the SneePath screening workflow, sending a text message with a link to download the WatchPAT mobile App.
- The progress bar updates according to each stage in the journey and to the office setting.
- You may convert a screening patient into a WP patient by either initializing a WP device or by changing the progress bar to Assign WatchPAT.

Screening Patient Card

3.5 Importing Patients

New patients can be added using the Import button which also allows creating multiple patients from a CSV file.

The Patient ID will be unique per Office. The following fields are supported:

- Patient ID (Mandatory)
- Office (Mandatory)
- Last Name
- First Name
- Referring Physician
- Gender
- Insurance Provider
- Other
- Status
- Custom Fields
- City
- ZIP
- Country

- Home Phone
- Country Code + Mobile Phone
- Email
- Number of Nights
- Bracelet Study
- Address
- State
- Work Phone
- Pacemaker
- Prefix
- DOB
- Logistics Comments
- STOP BANG
- Height
- Weight
- Epworth Score
- Neck Circ

3.6 First Time Installation

The first time you run the Initialization process or upload from the WP device, a one-time PC-App will be downloaded from CloudPAT.

CloudPATio (CloudPAT Input/Output) is a Windows application developed in Java and installed on the desktop computer for managing patient sleep data for Itamar Medical WatchPAT (WP) devices.

This includes the initialization and upload of studies and launch of the Study Viewer application.

Please note that the installation package size may reach a few hundred MBs since it includes some prerequisites such as Java JDK and AWS SDK (without installing them).

CloudPATio can be initially downloaded and installed as part of a required operation or deliberately through the Contact Us screen.

You should download the PC-App and run its installation - following the installation instructions.


Once it is installed – approve the use of the application as default.

CloudPATio provides two modes of installation:

1. A regular user installation making it operational only for the current Windows logged in user.
2. Admin user installation making it operational for any Windows logged in user. Admin installation shall require admin credentials for install, un-install and upgrade.


CloudPATio presents a taskbar icon when running allowing to terminate it from there.

Whenever there is a newer optional version, a suggestion to upgrade will be presented. Occasionally there may be a mandatory version that would have to be installed within a certain timeframe.




Download CloudPATio Screen

Double click the downloaded installation file.




Start the CloudPATio Installation


Follow the installation instructions:




Language Selection



Running the Installation Wizard



License Agreement Screen



Installation Completion Screen

3.7 Initialize WP200U/WP300 Device

Click on the Initialize Device button, select WatchPAT200U/300 option and input the desired settings on the screen:

- Bracelet Study (available only if office setting permits) - from this screen you may require a Tamper Proof Bracelet to be used as Positive Chain of Custody in Occupational Testing (e.g., truck drivers. See WatchPAT™ Operation Manual for further details).
- If the office setting permits a Multiple Nights Study: in order to record up to 3 nights with same WatchPAT device, select the number of nights needed. There are two options for Multiple Nights Study:
 - Two nights when no charging is needed between the nights (maximum length of study is approximately 8 hours)

- Three nights when charging is necessary between the nights (maximum length of study is approximately 10 hours)
- The Request Script option informs the Interpreting Physician that a script is requested.
- An associated interpreting physician may be selected for the study, and the option to change this selection will be allowed later.
- As part of the WP300 initialization process, there is an option to trigger a device test operation. This operation can also be invoked from the WP300 Operations Menu, by selecting the Device Test option.

Click on the Initialize Button once all inputs are set.

The system will prompt you to verify that the WatchPAT device is connected to the computer and then you should continue by selecting "Proceed with initialization".

A pop-up application will open and manage the initialization of the device.



Note

Make sure to connect the device, test the device and instruct the patient before giving out the WatchPAT™ device. Refer to the WatchPAT Operation Manual for more details.

When initializing a WP300 device, if an issue with the battery was detected, a warning message will appear. The user should consider putting in a new battery in this case.

The system will prompt you once the initialization process is completed.

If no Interpretation choice is listed (no physician selected) the report will not be sent for interpretation and a final report cannot be generated.

An office can configure a mandated interpretation selection as part of the device initialization process.

After the initialization is completed, a new entry is added to the Initializations section under the Studies & Initializations tab.

Overwrite protection: If a previous study was found on the WatchPAT™ device that was not uploaded to CloudPAT, you will not be permitted to overwrite that data. A prompt will appear: "Previous study found". In this case, you should upload the study before proceeding. This prevents you from inadvertently overwriting a study that was not uploaded.



Note

The office can update the patient's details as long as the report wasn't electronically signed. Once the report is electronically signed, the modified patient data will not be included in the report.

To prevent redundant work for the Interpreting Physician and the need to open a signed report, it is recommended not to assign the study until all patient and study fields are updated correctly.

Initialize Device

WatchPAT
ONE

☐

WatchPAT
300

☒

WatchPAT
400

☐

WatchPAT
ONE-M

☐

INTERPRETING PHYSICIAN

NUMBER OF NIGHTS

REQUEST SCRIPT
☐

BRACELET STUDY
☐

TEST DEVICE
☒

* Battery should remain in the device until study completion.

Cancel

Initialize

Initialize WatchPAT 200U/300 Device Screen


3.8 Initialize WatchPAT™ ONE Device


Click on the Initialize Device button, select WatchPAT ONE option and input the desired settings on the screen:


- Enter the WatchPAT™ ONE's Device serial number.
- Enter the PIN and PIN Type (hint) provided to the patient. Pin Type can be pre-set as an office setting.
- If there is a mobile number for the patient and the office setting for "Use last four digits of patient's mobile number" is enabled, the PIN would be auto populated from the last four digits of the patient's mobile number.
- Enter the max number of PIN Retries (default: 7, can be pre-set as an office setting).
- Accept default or change the Email for Notification, indicating the email for notification on study completion and upload to CloudPAT.
- Optional: Select the desired Interpreting Physician that the study will be sent to. If Interpreting Physician is not assigned, they can be selected after study is successfully uploaded to CloudPAT.
- An office can configure a mandated interpretation selection as part of the device's initialization process.
- Click on the "Initialize" button to complete the initialization process.


✕

Initialize Device


WatchPAT ONE


WatchPAT 300


WatchPAT 400


WatchPAT ONE-M

* SERIAL NUMBER

* PIN

**The last four digits of the Mobile Number

* MAX NUMBER OF PIN RETRIES

* PIN TYPE

Mobile Number ▼

**As configured for your office

* EMAIL FOR NOTIFICATION

Add

INTERPRETING PHYSICIAN

REQUEST SCRIPT

☐

Cancel

Initialize

Initialize WatchPAT ONE Screen

After the initialization is completed, a new entry is added to the Initializations section under the Studies & Initializations tab.

If patient has exceeded the maximum PIN retries, The counter may be reset by entering the Initialization section and selecting the “Reset Available PIN Retries” button.

If there is an error with the registration, you may select the Unregister button to cancel.

Initializations			
INITIALIZATION DATE	SERIAL #	DEVICE TYPE	STATUS
01-08-2021 12:44:04	XXXX	WatchPAT One	WP1 Registered

DEVICE TYPE

WatchPAT One

LOAD DATE

01-08-2021 12:44:04

Retest Device

PIN TYPE

Pin

PIN

1234

MAX NUMBER OF PIN RETRIES

3

Reset Retries

ASSIGNED PATIENT NAME

Mr. Erik Smith

Unregister

Reset Available PIN Retries

3.9 Initialize WatchPAT™ ONE Multi-Night Device

CP supports initializing multi-night WP1 devices (WP1-M) where the max Number of nights is 3 using a designated option in the Initialize Device screen.

The first Serial # will be the WP1-M, and the others will be of the additional probes.

- Optional: Select the desired Interpreting Physician that the study will be sent to. The study will be assigned to the interpreting physician as soon as the last study has been uploaded.

Once the Initialize button is clicked, the corresponding initializations shall take place and the PIN code shall remain the same for all initializations.

Initialize Device

WatchPAT
ONE

☐

WatchPAT
300

☐

WatchPAT
400

☐

WatchPAT
ONE-M

☒

* NUMBER OF NIGHTS

* PIN

**The last four digits of the Mobile Number

* MAX NUMBER OF PIN RETRIES

* EMAIL FOR NOTIFICATION

REQUEST SCRIPT

☐

* DEVICE SERIAL NUMBER

enter here the serial number that is presented on the box.

* PROBE2 SERIAL NUMBER

* PIN TYPE

**As configured for your office


Initialize WatchPAT ONE Multi-Night Screen

3.10 Initializing a WP400 device


Click on the Initialize Device button, select WatchPAT400 option and input the desired settings on the screen:

- Enter the WatchPAT™ 400's Device serial number.
- Enter Probe Serial Number.
- The PIN field is always the last 4 digits of the mobile number for WatchPAT400, and the PIN type is always Mobile Number. If there is no mobile number saved for the patient, in the WP400 Initialization screen a message will show saying that it is needed.
- Accept default or change the Email for Notification, indicating the email for notification on study completion and upload to CloudPAT.
- Optional: Select the desired Interpreting Physician that the study will be sent to. If Interpreting Physician is not assigned, they can be selected after study is successfully uploaded to CloudPAT.
- An office can configure a mandated interpretation selection as part of the device's initialization process.
- Click on the "Initialize" button to complete the initialization process.


Initialize Device




WatchPAT
ONE



WatchPAT
300



WatchPAT
400



WatchPAT
ONE-M

* SERIAL NUMBER

* PIN

**The last four digits of the Mobile Number

* EMAIL FOR NOTIFICATION

Add

REQUEST SCRIPT

☐

* PROBE SERIAL NUMBER

* PIN TYPE

Mobile Number

INTERPRETING PHYSICIAN

Cancel
Initialize

Initialize WatchPAT400 Screen

3.11 Device Initialization Details

In the Patient Card, under the “Devices and Initializations” tab, users can click on the caret for each initialized device to see the device details.

Here, information about the initialized devices can be seen, and additional functionality.

INITIALIZATIONS					
INITIALIZATION DATE	SERIAL #	PROBE SERIAL #	DEVICE TYPE	INTERPRETING PHYSICIAN	STATUS
▼ 05-15-2025 07:57:23	██████	██████	WP400		Registered
<div style="display: flex; justify-content: space-between;"> <div style="width: 30%;"> <p>DEVICE TYPE WP400</p> <p>LOAD DATE</p> <p>Retrieve Test</p> </div> <div style="width: 30%;"> <p>PIN TYPE Mobile Number</p> <p>PIN ██████</p> </div> <div style="width: 30%;"> <p>ASSIGNED PATIENT NAME</p> <p>Unregister</p> </div> </div>					
▼ 01-30-2025 08:11:16	██████		WP1		Completed
<div style="display: flex; justify-content: space-between;"> <div style="width: 30%;"> <p>DEVICE TYPE WP1</p> <p>LOAD DATE 01-30-2025 08:13:20</p> <p>Retrieve Test</p> </div> <div style="width: 30%;"> <p>PIN TYPE Mobile Number</p> <p>PIN ██████</p> <p>MAX NUMBER OF PIN RETRIES 7</p> <p>Reset Retries</p> </div> <div style="width: 30%;"> <p>ASSIGNED PATIENT NAME</p> <p>Unregister</p> </div> </div>					

Device Initialization Details Screen

3.12 Uploading a Study

This section covers uploading studies from the device or manually.



Note


WatchPAT ONE studies are automatically uploaded to CloudPAT once the sleep study is completed.

3.13 Normal Upload Mode (WP200U/WP300)

When a patient returns the WatchPAT™ device after a night study, simply connect WatchPAT™'s mini-USB port to any USB port on the PC by using WatchPAT™'s communication cable and click on the “Upload Study” button from any Screen.


Proceed by selecting the “Proceed with upload” button.

In case your browser asks for permission to run a command line for performing the upload, click on the button to allow – you may also select a checkbox to not receive this message again.



Allowing Windows Command Processor

During upload, a screen appears which details the various steps of the upload.



Patient Sleep Study Upload in Progress



Note

When a multi-night study is loaded all the night recordings are loaded automatically.



Note

Depending on your connection speed to upload, the number of nights included in the study and the type of device, the upload time will vary between 2 and 15 minutes.

Do not disconnect the WatchPAT device from the PC workstation until CloudPAT displays a message indicating that the study was successfully uploaded, otherwise the WatchPAT device might not be ready for initiation of a new WatchPAT study.

A notification message will be presented to the user as soon as the process completes.

✓ upload succeed! Redirecting to patient page...

Successful Upload Notification

Upon a successful upload, the system will present the study information under the patient's Studies & Initializations tab for observation of the study details.

Here you can set or change the desired interpreting decision or update the MIN DESAT FOR AHI/RDI if needed – both by clicking on the pencil icon for editing and on the checkmark for saving.

By the end of the WP300 initialization process, CP will notify the user in case there is an updated firmware version for this device.

By the end of the wp300 study upload process, if there is a newer firmware version, a suggestion will appear to perform the upgrade straight from CloudPAT. Once accepted, it will trigger the upgrade operation. This operation can also be invoked from the WP300 Operations using the Firmware Upgrade option.

If you attempt to upload a WP300 study to a patient who already has a study with the same date, a message will be displayed.

STUDY DATE	SERIAL #	AHI	ODI	VALID SLEEP TIME	# OF NIGHTS	ANALYSIS READY	INTERPRETING PHYSICIAN	MIN DESAT FOR AHI	MIN DESAT FOR ODI	REPORT READY
> 08-11-2022	[REDACTED]			0 hrs, 0 min	1	Yes		3%	4%	No

INITIALIZATION DATE	SERIAL #	DEVICE TYPE	INTERPRETING PHYSICIAN	STATUS
08-11-2022 15:51:34	[REDACTED]	WatchPAT 300		Log Edit
08-11-2022 15:48:52	[REDACTED]	WatchPAT 300		

In-Row Editing in Study Details Screen

3.14 Manual Mode (WP200U/WP300)

In the rare event that a “manual” upload is needed (such as an internet outage and the sleep study was uploaded locally using zzzPAT), select “Switch to Manual Upload”.

The study files must be zipped together into one zip file to be selected through a file browser.

These files are patient.dat, sleep.dat and if existing, also log.dat, sleep1.dat and sleep2.dat (for multi night studies).

Click on the Upload Area for opening a file browser to locate the desired zip file.

Once added, click on the Upload ZIP File button to continue.

From this point on, the process is similar to the Normal Upload Mode described above.

On attempt to load a study that has been initialized not in CP 3.0, while the user is connected to more than one office, CP will present a dropdown of the user's connected offices, allowing the user to select the desired office and then create the patient and the study under this office.

Upload Study Screen

3.15 Manual Upload of WatchPAT ONE Study

In the rare occasion, the WatchPAT ONE (WP1) study was not automatically uploaded to CloudPAT after the sleep study recording was completed, you may upload the study manually via the Initializations section.

Select the relevant WP1 registration entry and click on “Retrieve Test” button.

You may continue working with CloudPAT and the study will be uploaded in the background.

An email will be sent to the email address provided during the WP1 registration when the study is successfully uploaded to CloudPAT.

Manual Retrieval of WP1 Study

3.16 Searching Studies and Reports

Selecting the “Patients” icon opens a search screen where you can view all or a subset of the studies and reports available to the current user. Clicking on the table's icon will perform a table refresh.


The following filters are available: Patient ID, First Name, Last Name, Email, Office, Status, Custom Status, Period, Sleep Apnea, Serial #, Study Date, Study Load Date, Study Close Date, Waiting for Interpretation, Interpreting Physician, Patient Created Date, Referring Physician, Insurance Provider, and Report Viewed (indicates if the report file has been opened by a provider from the office).

When the filter fields include study information, the result is a list of patients having at least one study matching the search filter.

The table including the search results is customizable and includes the majority of the search fields.

Search results can be exported to CSV files.

- Export to CSV: Including patients' details.
- Export Summary Report: Including some analysis information.
- Export History: Including progress bar changes.



Viewing Sleep Study Reports

3.17 Study Management

Once a study is uploaded and until the Interpreting Physician locks the report, there are several tasks that a user can do.

3.17.1 Understanding the Study Screen

This screen presents information about the study.

The more frequently used fields are displayed on the collapsed row.

Changing the Minimum % of Desaturation for AHI/RDI Calculation and Assigning Interpreting Physician can be done here, using the pencil icon for editing, and the check mark icon for saving.

Selecting the caret expands the view to include all additional information from the study: errors or warnings from the analysis, Bracelet information, and if the Request Script box is checked or not.

If the office settings allow it, you can see the Analysis File under “Files”.

The office settings control what information in the study is shown in the row.



Note

To prevent redundant work for the Interpreting Physician and the need to open signed reports, it is recommended not to assign the study until all patient and study fields are updated correctly and all the relevant files are attached.

The screenshot shows the CloudPAT Patient Card interface. At the top, there's a progress bar with five steps: Created, In Testing (initiated), Study Loaded, Sent to Interpretation, and Study Completed (5). Below this, patient information is displayed in a form: PATIENT ID (ronen111), NAME (Ronen), OFFICE (Ronen Office), GENDER, DATE OF BIRTH, AGE, REFERRING PHYSICIAN, COUNTRY CODE, MOBILE PHONE, and EMAIL. There are tabs for Studies & Initializations, Details, Diseases & Medications, Attached Files, and App Questionnaires. The 'STUDIES' section shows a table with columns: STUDY DATE, SERIAL #, AHI, ODI, VALID SLEEP TIME, # OF NIGHTS, ANALYSIS READY, INTERPRETING PHYSICIAN, MIN DESAT FOR AHI, MIN DESAT FOR ODI, and REPORT READY. The first row shows data for 08-14-2022. Below the table, there are fields for ASSIGNMENT DATE, LOAD DATE, LOCK DATE, BRACELET STUDY, BRACELET ID, REQUEST SCRIPT, and ERROR / WARNINGS. At the bottom, there are buttons for RAW Data, RTF, Log, and Edit.


STUDY DATE	SERIAL #	AHI	ODI	VALID SLEEP TIME	# OF NIGHTS	ANALYSIS READY	INTERPRETING PHYSICIAN	MIN DESAT FOR AHI	MIN DESAT FOR ODI	REPORT READY
08-14-2022	[REDACTED]	22.4 (3%)	15.0 (4%)	5 hrs, 19 min	1	Yes	sp Ronen	3%	4%	Yes

Study Details Screen

3.17.2 Attaching Files to a Patient

Once a patient is created, the user can choose to add documents to the patient record. These documents will be available for the Interpreting Physician to review when preparing the report.

To attach a file – enter a description and select the file to be uploaded using the upload area.




Attach File Screen

The user can also remove or download the attached files by clicking on the respective buttons.

Offices can mandate to configure at least one attached file before the initialization process.

Additions, updates or deletions of attached files are tracked in the Audit Trail Screen which may be viewed by clicking on the Log button.



Deleting Attached Files

3.17.3 Changing Interpreting Physician

The user can assign a non-assigned study to any of the Interpreting Physicians connected to the office. The user can also change the assignment between the different Interpreting Physicians connected to the office.

Assigning interpreting physicians can be done through the study's row using the pencil icon for edit and check mark icon for saving.

Notification of a new study will be sent to the assigned Interpreting Physician while another notification will be sent to the previous Interpreting Physician indicating that the study has been reassigned.

Offices may turn on an Eligibility Check for Interpreting Physicians. This would compare the states where the interpreting physician is board certified in with the Patient's state of residence and would allow a study assignment only if they match.

3.17.4 Interpretation Pools

Interpretation Pool is a group of interpreting physicians that can be assigned to an office; studies may then be assigned to a pool instead of a designated interpreting physician.

Each pool member could then pull a study from the pool and become the designated interpreting physician for this study.

Assignment to an interpretation pool is done similarly as with interpreting physicians.

An interpreting pool member may return a study to the pool for someone else to proceed with, as long as the study isn't locked by clicking on the release button.

A pool member may be configured as a pool manager and as such they would be able to re-assign studies and get daily emails about studies that have not been assigned.

A pool member may be configured as not authorized to lock a study, only to save changes to it.

If the Eligibility check is configured for the office, an interpretation pool could be assigned only if it includes at least one member that is certified to interpret in the patient's state and only certified members could pull the study.

Whenever a pool member is configured as the Default Assignee, they would be assigned to any new study assigned for this pool, getting these studies added to their Studies table. The default assignee shall be able to return studies back to the pool for other members to pull them.

3.17.5 View Analysis File

Some offices have settings that allow their users to view the analysis before the report is ready.

This setting is usually given to sleep labs using the CloudPAT where the expertise needed to understand the results are part of the office skill set.

Users of such offices can click the analysis file link immediately after uploading the study.

When a multi-night study is uploaded, all the results from the reports will appear in an analysis file, one after another, followed by a summary page.

Offices have the option to enable a setting to allow providers to download the statistics.XML file (including the full set of the analysis information).

Offices also have the option to enable downloading EDF file of an analysis.

**Note**

An analysis file is created during the analysis and may include patient details entered at the time of the analysis. Patient details will not change in the report until the report is locked or re-run.

3.17.6 Request for a Script

In case the Interpreting Physician is required to add a script to the report, the “Request Script” box should be checked. Save the study by selecting the “Save” button, and complete the request by sending it to the Interpreting Physician.

The Interpreting Physician will get a reminder to add a script upon locking the report.

3.17.7 Changes of Study Parameters

The user can change patient information and some study parameters as long as the report is not closed.

Some of the parameters will require running the analysis again:

- AHI index – according to 3% or 4% desaturations.
- Pacemaker.
- Gender
- Date of Birth

After making changes to these parameters, the user should rerun the analysis to get the new results. A reminder about the need to rerun the analysis will appear on the page together with an option to rerun the analysis. A re-analyze notification will be sent to the Interpreting Physician.



Note

Per office setup, Desaturation Threshold may automatically be set to 4% for patients insured by Medicare.

A warning message may appear if the Desaturation Threshold is not 4% for Medicare insured patients.

The screenshot shows the CloudPAT Patient Card interface. At the top, there's a progress bar with five steps: 1. Created (08-11-2022), 2. In Testing (initiated) (08-11-2022), 3. Study Loaded (08-11-2022), 4. Sent to Interpretation, and 5. Study Completed. Below the progress bar, patient details are listed: Patient ID (ronen special ch), Name (ronen special ch), Office (Ronen Office), Gender, Date of Birth, Age, Referring Physician, Country Code, Mobile Phone, and Email. A warning message in an orange box states: "Rerun is required. Analysis is outdated as: The percentage of Min Desat AHI and/or ODI has been changed." Below this, there are tabs for Studies & Initializations, Details, Diseases & Medications, Attached Files, and App Questionnaires. The 'STUDIES' tab is active, showing a table with columns: STUDY DATE, SERIAL #, AHI, ODI, VALID SLEEP TIME, # OF NIGHTS, ANALYSIS READY, INTERPRETING PHYSICIAN, MIN DESAT FOR AHI, MIN DESAT FOR ODI, and REPORT READY. The first row shows a study from 08-11-2022 with a serial number of 105904, 0 hrs, 0 min valid sleep time, 1 night, and analysis ready status. Below the table, there are fields for ASSIGNMENT DATE, LOAD DATE (08-11-2022), and LOCK DATE. There are also checkboxes for BRACELET STUDY, BRACELET ID, and REQUEST SCRIPT. An 'Inconclusive REM. Incomplete analysis due to a Noisy PAT signal.' message is displayed. At the bottom right, there's a 'Log' button and an 'Edit' button. A warning message at the bottom right states: "Rerun is required. Analysis is outdated as: The percentage of Min Desat AHI and/or ODI has been changed."

Rerun Analysis Indication

3.18 Using CloudPAT as an Office Admin

The Office Admin is a role designed to manage and administrate the Office, and has further capabilities without relying on Itamar Admin for operations.

An Office Admin has the following capabilities and permissions:

Creating office Providers:


An Office Admin can create Office Providers by navigating to the "Access" tab and clicking the "New" button, then fill the Provider Details fields.

Creating Office Interpreting Physicians:

An Office Admin can create Office Interpreting Physicians by navigating to the “Access” tab and clicking the “New” button, then fill the Interpreting Physician Details fields.

Creating Office Sleep Physicians:

An Office Admin can create Office Sleep Physicians by navigating to the “Access” tab and clicking the “New” button, then fill the Sleep Physician Details fields.



“Create New” dropdown Menu

Office admins have full control over SleepPATH settings.

3.19 SleepPATH Questionnaires

App Questionnaires are introduced to upgrade the HSAT experience. On a digital version, these questionnaires are customizable by provider and patient and help sleep physicians quickly and efficiently centralize data management and enhance the diagnosis process.

The SleepPATH settings can be configured by an Itamar administrator or by an Office Admin user granted with this permission, also including the possibility to create Y/N custom questions for both screening and the study questionnaires, to be added in any of the sections. Changes in the SleepPATH settings are captured in Audit Trail screen.

The provider can turn questionnaires On/Off from the App Questionnaire tab and select elements from Bedtime, Morning, Assessment Questions, Risk Assessment (STOP BANG) that will later show in the patient’s mobile application.


Once questionnaires are turned on, the questionnaire menu is displayed.

When using questionnaires, Email and Mobile fields become mandatory and CP issues an error message on attempt to Initialize a Device if empty.

The Patients Table includes an icon to indicate the App Questionnaires status (In Progress, Completed, Aborted, NA).

When selecting the "Resend Entry Code" button an SMS message will be sent to the patient's mobile phone with a new entry code to be used in the APP.

The answers are added to the sleep report and are also available in the Studies & Initializations tab.



App Questionnaires

Office Admins have full control over sleepPATH settings through the "sleepPATH" Tab under Office settings.

3.19.1 Learning center:

Learning center is a mobile APP service that allows the office to send the patients educational information pre and post sleep study and to get confirmation that the patient viewed the information.

Pre-Study:

The pre-study education process is triggered during the device initialization process.


It may include the following items: Educational Videos, About sleep apnea, Common symptoms, Associated risks, Treatment options, FAQs, Healthy sleep habits.

Once checked by the provider, the relevant content will be available for the patient in the mobile application.

Post-Study:

The post-study education process may be triggered or the interpreting physician after the study has been locked in a similar fashion.

It may include the following items: Educational Videos, Sleep Apnea, Common Symptoms, Associated Risks, Treatment Options, FAQs, Patient Preferred Therapy (Survey), Healthy sleep habits and the patient's treatment preference.



SleepPATH Learning Center

3.20 IDTF Customers - HSAT Order Form

There exists an optional tab called “HSAT Order Form” if the office is a ZOLL IDTF customer submitting orders using CloudPAT. This tab enables the creation of a sleep test order in a PDF format.

Any of the Patient Information fields that were entered/modified in the e-Order tab will update the Patient Details tab and vice versa.

There are several sections, and each section holds different fields. Certain logic is applied to the form to ensure the IDTF has the information it needs to process the order.

The logic for submitting an E-order: Once the E-Order button is pressed, CP executes a validation of input data, then the user shall put in the “Prescriber Login” field the User Name of the logged in-user, and in the “Prescriber Password” field the password of the logged in user. A PDF is then created with the entered data under the Attached Files tab.

You may attach PDF files to be included in the output PDF.

Orders can be updated and repeated.

If the “Path to Therapy” setting is turned on, the form will include a checkbox called “Information Sharing”. If turned on, a selection of a therapy provider from the dropdown will be mandatory for producing the pdf. This will send the study to the chosen therapy provider.

The screenshot displays the CloudPAT Patient Card interface. The top navigation bar includes 'CloudPAT', 'Patients \ Patient Card', and a 'WatchPAT Arrhythmia Detection' banner. The main content area shows the 'HSAT Medical Order Form' with tabs for 'Studies & Initializations', 'Details', 'IDTF', 'Diseases & Medications', 'Attached Files', 'SleepPath™', and 'HSAT Medical Order Form'. The form includes a table with columns 'FILE NAME', 'UPLOAD DATE', 'UPLOADED BY', and 'CREATION REASON'. Below the table is a 'Create Order Form' button and a 'Creation Reason' dropdown menu. A note states: '**Please make sure all mandatory fields are filled, if some are disabled please turn to your system admin and make sure he updates the missing fields**'. The 'PRESCRIBER INFORMATION' section contains fields for 'Prescriber Name', 'NPI', 'Prescribing Facility', 'Address', 'City', 'State', and 'Zip'. At the bottom right are 'Cancel', 'Save & Initialize', and 'Save' buttons.

HSAT Order Form

3.21 Viewing Reports

After a patient sleep study has been sent to your Interpreting Physician for review and interpretation, and the Interpreting Physician has reviewed and finalized the report, the user who uploaded the study will receive an email notification that your patient sleep study report is available.

Alternatively, you may poll the system by pressing on the caret of the study and checking the “Report Ready” field.

If yes is shown, you can review the signed sleep study report and the attached files by clicking on the PDF File button.

The report will include comments and diagnosis from the Interpreting Physician. PDF file attachments added by the Interpreting Physician will be appended at the end of the report. Other attachments (non-PDF) will only be viewable by downloading them (clicking on the link to the file).

In case of a multi-night study, the analysis of each night will appear one after another in the report, followed by a summary page at the end.

Subject to office configuration, the report may show the other AHI / ODI values in the report based on the alternative MIN DESAT values.

Based on the office setting, it is possible to get the report in a secondary format which is either HTML or RTF.

If configured for the office, the report will include a page with OXI and Pulse Rate Histograms.

If configured for the office, the report will include Arrhythmia events.

The screenshot displays the CloudPAT Patient Card interface. At the top, a progress bar shows the study status: Created (08-11-2022), In Testing (initiated) (08-11-2022), **Study Loaded** (08-11-2022), Sent to Interpretation, and Study Completed. Below this, patient details are listed: Patient ID (ronen special ch), Name (0009kkceee660e), Office (Ronen Office), Gender, Date of Birth, Age, Referring Physician, Country Code, Mobile Phone, and Email. A sidebar on the left contains navigation icons. Below the patient details, there are tabs for 'Studies & Initializations', 'Details', 'Diseases & Medications', 'Attached Files', and 'App Questionnaires'. The 'Studies & Initializations' tab is active, showing a table of studies. The table has columns for Study Date, Serial #, AHI, ODI, Valid Sleep Time, # of Nights, Analysis Ready, Interpreting Physician, Min Desat for AHI, Min Desat for ODI, and Report Ready. One study is listed for 08-11-2022 with a serial number of 0000, 0 hrs, 0 min valid sleep time, 1 night, and analysis ready status. Below the table, there are fields for Assignment Date, Load Date (08-11-2022), and Lock Date. To the right, there are checkboxes for 'BRACELET STUDY' and 'REQUEST SCRIPT'. An 'ERROR / WARNINGS' section shows a red message: 'Inconclusive REM. Incomplete analysis due to a Noisy PAT signal.' At the bottom, there are buttons for 'RAW Data', 'RTF', 'Log', and 'Edit'.

Download Buttons for the Reports



Note

Contact Itamar Medical Support for the availability of other file type options for reports.

3.22 Managing Devices

This is an inventory management feature allowing the users to keep track on the multi-usage devices of the office.

A device can be added through the new button, or once it is initialized through the Initialize Device button.

3.22.1 Devices Table

You can search devices from the table using designated filter fields:

Serial #, Device Type, Status, Office, Status Change Date, Device Assignee, Firmware Version.

Once a study is being uploaded from the device, its status changes to Available.

SERIAL #	DEVICE TYPE	OFFICE	DEVICE ASSIGNEE	STATUS	STATUS CHANGE DATE	FIRMWARE VERSION
[REDACTED]	WatchPAT 300	Ronen Office	[REDACTED]	Available	08-23-2022 18:12:08	
[REDACTED]	WatchPAT 300	Ronen Office	[REDACTED]	Assigned	08-17-2023 14:39:46	4.2.1023
[REDACTED]	WatchPAT 300	korenNEW	[REDACTED]	Assigned	10-31-2023 13:58:29	
[REDACTED]	WatchPAT 300	korenNEW	[REDACTED]	Out of Service		
[REDACTED]	WatchPAT 300	korenNEW	[REDACTED]	Assigned	08-22-2022 04:15:30	
[REDACTED]	WatchPAT 300	korenNEW	[REDACTED]	Assigned	08-22-2022 04:35:43	
[REDACTED]	WatchPAT 300	korenNEW	[REDACTED]	Assigned	08-22-2022 04:11:07	
[REDACTED]	WatchPAT 300	korenNEW	[REDACTED]	Available		
[REDACTED]	WatchPAT 300	korenNEW	[REDACTED]	Assigned	07-20-2022 21:11:16	
[REDACTED]	WatchPAT 300	korenNEW	[REDACTED]	Assigned	07-20-2022 21:36:47	
[REDACTED]	WatchPAT 300	korenNEW	[REDACTED]	Assigned	07-20-2022 23:25:02	
[REDACTED]	WatchPAT 300	korenNEW	[REDACTED]	Assigned	08-22-2022 03:41:32	
[REDACTED]	WatchPAT 300	korenNEW	[REDACTED]	Assigned	08-22-2022 06:51:44	

Devices Table Screen

3.22.2 Device Details

Selecting the Device row, opens the Device Details Screen showing additional information.

The screen includes the following fields: Device Type, Status, Comment, Bracelet Study, Serial #, Assigned Patient ID, Office, Firmware Version, History.

A device may be labeled as “Out of Order” by selecting the “Out of Order” button.

Initialization of a device that is marked as “Out of Order” is prohibited.

A device can be transferred between offices by using the Office dropdown.

The screenshot shows the 'Device Details' screen with the following fields and values:

- DEVICE TYPE ***: WatchPAT 300 (with a device icon)
- STATUS**: Available
- COMMENTS**: (Empty text area)
- BRACELET STUDY**: No
- SERIAL # ***: 105904
- ASSIGNED PATIENT ID**: No
- OFFICE**: (Dropdown menu showing 'Ramat Hashikma')
- FIRMWARE VERSION**: 4.2.1023
- HISTORY**: A list of three entries, each showing a date (09-03-2023), a status (Study), and a time range (10:00-23:02-23).

At the bottom, there are three buttons: a pink button with a 'G' icon, a grey 'Cancel' button, and a green 'Save' button.

Devices Table Screen

3.23 Managing Referring Physician

It is possible to manage referring physicians at the office level without waiting for system administrators.


A referring physician may be added, edited, or named as inactive by the provider.

3.23.2 Referring Physician Details

Selecting the Referring Physician row opens the Referring Physician Details Screen displayed with additional information.

A Referring Physician can be designated as Inactive by using the designated checkbox.

If the user is associated with multiple offices, they may move a referring physician to one of the other offices.



The form is titled "Referring Physician" and contains the following fields and controls:

- PREFIX:** A dropdown menu with "Dr." selected.
- FIRST NAME:** A text input field containing "Sam".
- MIDDLE INITIAL:** A text input field containing "J".
- * LAST NAME:** A text input field containing "Shoval".
- SUFFIX:** A dropdown menu with "MD" selected.
- JOB TITLE:** An empty text input field.
- * OFFICE:** A dropdown menu with "office" selected.
- EMAIL:** An empty text input field.
- * ACTIVE:** A toggle switch that is currently turned on (blue).
- SPECIALTY:** A dropdown menu with "Cardiologist" selected.
- NPI / LICENSE:** A text input field containing "1234567890123".
- INSTITUTE:** An empty text input field.
- ADDITIONAL PHONE NUMBER:** An empty text input field.
- MOBILE PHONE:** An empty text input field.
- Buttons:** "Cancel" (grey) and "Save" (green) buttons at the bottom.

Referring Physician Details Screen

3.24 Managing Lists

Managed Lists include medications, diseases, statuses, and insurances.

The user can add items to the list and edit or delete them.

Selecting the "Managed Lists" icon opens a search screen where you may search for items globally or locally. The Global list cannot be changed.

It is possible to create multiple items using the semicolon separator.

NAME	OFFICE	TYPE
United Health		INSURANCE
Thyroid hormone replacement		MEDICATION
Statins		MEDICATION
Stable Angina		DISEASE
Snoring		DISEASE
Sleep apnea		DISEASE
Schizophrenia		DISEASE
romen status	office	CUSTOM_STATUS
romen ins	office	INSURANCE
RLS		DISEASE
PLMD		DISEASE
Peptic ulcer		DISEASE
Patient Pay		INSURANCE
Oxygen therapy		MEDICATION
Other Private		INSURANCE
Osteoporosis		DISEASE
office_1_medication1	office	MEDICATION
office_1_disease1	office	DISEASE

Managed Lists

3.25 SleepPATH Communications Tab

If a patient belongs to an office which has enabled "Active Communications", an additional tab will be added to the patient card, called "Communication", under "SleepPATH" tab.

This service sends texts to the patient, which are daily reminders, reminding the patient to complete the study.

When Custom Status is set to 'Patient Reject Test' or when the office is connected to PCP Referral Sales Region, no further communication shall be performed for that patient.

CloudPAT Patients \ Patient Card

Studies & Initializations | Details | WPT | Diseases & Medications | Attached Files | **SleepPATH™**

Questionnaires | Learning Center | **Communication**

SleepPATH™ \ Communication

✓ REMINDER MESSAGES (Automated)

- These are automated messages sent to patients that have been assigned a WP and have not yet completed the study.
- These messages will be sent until they are disabled or reach the expiration date.
- Patients cannot reply to these messages.

Sending Frequency Every Days
(Minimum every 5 days, Maximum every 10 days)

Expiration Last Sending Date - After Times
(Maximum after 10 times, Minimum after 2 times)

Text - Complete the WP1 Study reminder

On behalf of korenNEW, we remind you to complete your Sleep Study. Please click the link below to START the process: \$(link). Thank you.

Sent: -
Next: -

Patient Communication Tab

3.26 Using CloudPAT as a Therapy Provider

The Therapy Provider is a role designed for sleep physicians for treating patients.

Users who are defined as Therapy Providers will have an additional menu item called “Therapies” which leads to the therapies table.

The therapies table contains all the patients that were assigned to the therapy provider, and their studies have been locked.

There, the therapy provider can choose to accept or reject the therapy. This is done by clicking the “Accept Therapy” or “Decline Therapy” buttons respectively.

Clicking on the “Accept Therapy” button, then clicking on the patient’s row again will navigate into the therapy screen.

Clicking on the “Decline Therapy” button will remove the therapy from the therapy provider’s pending therapies.

ACCEPT THERAPY	PATIENT ID	OFFICE	REFERRING PHYSICIAN	INTERPRETING PHYSICIAN NAME	STUDY DATE	STUDY LOAD DATE	VALID SLEEP TIME	RDI
✓ Accept Therapy	[REDACTED]	[REDACTED]	[REDACTED]	[REDACTED]	05-06-2025	05-06-2025	5 hrs, 46 min	36.9
✓ Accept Therapy	[REDACTED]	[REDACTED]	[REDACTED]	[REDACTED]	05-06-2025	05-06-2025	5 hrs, 46 min	36.9

1
Total 2

Therapies Table

4 Using CloudPAT as an Interpreting Physician

4.1 Main Access Screen

When you sign in CloudPAT as an Interpreting Physician you will be presented with the Interpreting Physician's Home screen.

Here are the main actions that can be selected from this screen:

1. Review studies – reviews new or display completed studies.
2. Manage templates and NPIs – defines new or edits existing NPIs and templates used for Clinical diagnosis and recommendations.
3. Update settings including Interpreter Signature.

PATIENT ID	PATIENT	OFFICE	REFERRING PHYSICIAN	STUDY DATE	STUDY LOAD DATE	VALID SLEEP TIME	AHI	ODI	WAITING FOR INTERPRETATION
ronen100		office		04-04-2022	04-04-2022	5 hrs, 19 min	22.40 (3%)	15.00 (4%)	
1234567890123456	1234567890123456, 1234567890123456	office		04-21-2021	04-21-2021	4 hrs, 58 min	24.10 (3%)	16.10 (4%)	482
ronen001		office		01-02-2020	04-11-2021	4 hrs, 58 min	23.90 (3%)	16.10 (4%)	
ronen 110		office		02-09-2021	02-09-2021	0 hrs, 0 min			
cp804		Ronen Office		01-02-2020	01-28-2021	5 hrs, 22 min	22.20 (4%)	15.10 (4%)	138
cloud_4802		Ronen Office		01-25-2021	01-25-2021	0 hrs, 0 min			567
cloud_4801	David	Ronen Office		01-02-2020	01-25-2021	4 hrs, 58 min	23.90 (3%)	16.10 (4%)	

Total 7

Interpreting Physician's Home Screen

4.1.1 Studies Assigned to the Interpreter

After a patient's sleep study has been sent for interpretation, you will receive an email notification that a sleep study is awaiting interpretation. Alternatively, you may poll this by selecting "Unlocked" Status in the filter of the main screen which will display all unlocked studies waiting for you. You may also select the "Locked" status for viewing completed (signed) sleep study reports. You have the option to review subsets of your studies by selecting a particular Office, Patient name, etc.

CloudPAT supports an option to send Interpreting physicians a summary of their pending studies through a daily email.

**Note**

Click on the CSV icon if you want to obtain an Excel file containing a tabular presentation of the search results.

4.1.2 Reviewing a Study

When you filter for “Unlocked” status, select the study you want, and the screen will display the study details for you to review. Here are the main actions that can be selected from this screen:

1. Interpretation Details tab including:

- Input fields for the interpreter to fill in the Diagnosis and Recommendations. You may use the Diagnosis and Recommendations templates if you wish to insert an already pre-defined paragraph in the corresponding section. Multiple templates can be used.
- Errors/Warnings – in case of errors or warnings returned from the analysis, a message will appear in the Errors/Warnings box.
- Change Min Desat for AHI/RDI – change the value used for analysis. Change becomes effective only after saving. If this occurs, a message will display indicating that rerun analysis is required for the new setting. Click on “Rerun Analysis” to see the results from the new threshold.
- NPI – you can add the NPI to the report if have one. If you have more than one NPI, select the right NPI number that matches the office details.
- Add attachments like prescriptions or any other files to the report. If you add PDF files, they will be added to the report itself. Any other files will be available to download in the office. To attach a file – enter a description and select the file to be uploaded using the upload area. The user can also remove or download the attached files by clicking on the respective buttons.
- Interpretation History section listing interpretation related operations done on the study
- Interpretation Comments section facilitates communication between interpreting physicians and pool members about the study.
- Interpreting Physicians and Sleep Physicians can delete their Interpretation Comments.
- A Sleep Physician can navigate from the Study Details screen back to the respective Patient Details screen by using a button called “Go to the Patient Record”.
- It is possible for the interpreting physician to add a statement to the sleep report that the study was manually scored / analyzed alongside with the interpreting physician name by clicking on a designated button.

2. Therapy tab:

- Allows the interpreter to produce a Therapy Prescription which when enabled retrieves patient information which can be edited and prescription items selected from the list.
- Once the Generate Script button is clicked, an auto-generated PDF file shall be created and displayed, including the submitted information to be attached to the report during the Save and Lock operation.
- The supported therapies are CPAP and Oral Appliance.

CloudPAT®
By Itamar

Studies \ Patient

Contact Us

Run Analysis

Preview

Open in Viewer

Download Study File

Interpretation Details

Therapy

Patient Details

Office Details

Interpreter Attached Files

Enable Script Prescription

Create Script

Patient Information

FIRST NAME*

LAST NAME*

PATIENT ID

GENDER

DATE OF BIRTH

BMI

STOP BANG SCORE

PATIENT MOBILE PHONE

NECK CIRC. (CM)

Prescription

DIAGNOSIS*

LENGTH OF PRESCRIPTION*

TREATMENT TYPE*

PATIENT EDUCATION

CLINICAL NOTES

HUMIDIFIER*

CPAP MASK / INTERFACE / DELIVERY SYSTEM*

Interpreting Physician Details

PHYSICIAN NAME*

NPI*

LICENSE

ADDRESS

STATE

CITY

ZIP CODE

PHONE

FAX

PRODUCE PRESCRIPTION

PRODUCE SCRIPT

Cancel

Save

Save & Lock

Therapy Tab

3. Patient Details Tab– Patient details are displayed.
4. Office Details Tab – by clicking the “Office Details” button, you can see the details of the office that sent the sleep study: address, phone and minimum % Desaturation for AHI/RDI Calculation.
5. Interpreter Attached Files Tab: Showing files that were attached by the interpreter.
6. Diseases & Medications tab: Showing Diseases, Medications and Clinical Comments entered by the provider.
7. Preview Report Button- You can click the “Preview Report” button to see the report that was produced automatically by CloudPAT. According to the Respiratory Indices, Various Statistics, and Sleep Stages, you will make your interpretation and fill in the Diagnosis and Recommendations sections. This also includes the patient information and your last saved interpretation.
8. Open in Viewer button (RECOMMENDED) – This utility avoids the need to install zzzPAT on the PC, to download the study to the PC and upload the results. If you wish to see more details on the study’s signals or to review and edit the events, you can use the “Open in Viewer” function to start the zzzPAT utility viewer on your local computer. Refer to the zzzPAT Operation Manual for a description of the operational capabilities of the Viewer.



Note

When you first use “Open in Viewer” on your PC, installation of a small PC-app, CloudPATio will begin (if not installed already as part of the device initialization / study upload process).

9. Download Study Files Button (Legacy, only if enabled to the office) – Manually download the study to zzzPAT, and manually upload editing results to CloudPAT. This legacy option is kept in CloudPAT during the migration to Open in Viewer method.



Note

Make sure to upload the correct manual event file. The filename will be comprised of

ID_[patient ID]_

SD_[study Date & time]_

events_analysis_user.zpt

other file will generate an error during the post upload analysis and will not change the analysis file

10. Save – saves your input without locking the report.
11. Save and Lock Report – When you are satisfied with your interpretation select the “Save and Lock Report” button. At this point the study report will be automatically digitally signed with your name, tagged as completed and be available for the Physician’s review.
12. In order to comply with part 11 regulations, an office may require a re-entry of the password upon locking the study.



Note

Report can be locked only if created during updated analysis. If analysis is outdated, the system will display a message and the “Save and Lock Report” will be disabled. After rerunning the analysis the “Save and Lock Report” button will be available again.

13. Unlock a Report – interpreting physicians can unlock their own locked studies by clicking on the Lock icon on the List of Studies Screen. An Audit Trail Record will be added to the study indicating the user and when the previously locked report will be added.
14. Notifications – you may receive notifications which will be displayed on the screen in the following cases:
 - If a script was requested by the office, a “Study requires a script” notification will be displayed on the screen. Also, a reminder will appear before locking the report.
 - If the office changed one of the following analysis parameters: Minimum % Desaturation for AHI/RDI calculation, Pacemaker, Gender, Date of Birth, or a rerun of the analysis is required to update the analysis result. In case rerun was not initiated by the office, a notification that rerun is required will be displayed on the screen and the “Rerun Analysis” button will appear. You won’t be able to lock the report before re-running the analysis.
15. Summarized PDF report – CloudPAT enables to produce a patient letter based on the PDF report through the WP Interface. This enables a button in Studies and Initializations called “Patient Letter”, and clicking on it will open a more summarized version of the PDF report.
this feature is enabled only for offices that turned it on through Office Settings.
16. Re-generate report after patient changes - CloudPAT enables to re-generate the Sleep Report through the WP Interface in cases where edits were made to the patient details (e.g., height, weight, age, name, DOB etc.) with no need to re-run the analysis.

Interpretation Details Screen

4.1.3 Viewing a Completed Study

When you filter for “Locked” status, select the study you want, and the study details screen will appear.

Here are the main actions that can be selected from this screen:

1. View report - you can click on the “View report” button in order to see the full final report (with your interpretation included).
2. Open in Viewer (read only).
3. View data on all tabs.

A Completed Study Screen

4.2 Manage Templates & NPIs

This section explains how to manage the lists of templates and NPIs available to each Interpreting Physician. Selecting the “Settings” icon will open a table used to browse for templates and NPIs. You can narrow the search by using the “Type” field.

NAME	TYPE	TEXT	ACTIVE
PHB001	NPI		Yes
Reco	RECOMMENDATION	Reco	Yes
CPAP Recommendation	RECOMMENDATION	If a patient is not able or willing to undergo CBT-I, then short-term (3-6 months) pharmacotherapy, complementary and integrative health (CIH) treatments, or both can be considered. If pharmacotherapy is considered, a thorough sleep history and evaluation for contraindications for pharmacotherapy are recommended. Nonpharmacological interventions are also encouraged. If pharmacotherapy is initiated, low-dose doxapram or nonbenzodiazepine benzodiazepine-receptor agonists (NBRAs) can be considered, though the excess risk for complex sleep-related behaviors with the NBRAs should also be discussed.	Yes
depression	DIAGNOSIS	As a medical diagnosis, depression denotes a prolonged state of dejection, emotional emptiness, self-recrimination, lack of energy, apathy, anxiety, sleep disorders, and lack of appetite.	Yes
Get CPAP	RECOMMENDATION	The following is the CPAP script:	Yes
90210	NPI		Yes
ronen diag	DIAGNOSIS	ronen diag details	Yes

Templates & NPIs Screen

4.2.1 Templates

Templates are predefined paragraphs ready for use while reviewing a study and entering the Clinical Diagnosis or Recommendations sections. For example, you may define a diagnosis template for Severe Sleep Apnea and a different one for Mild Sleep Apnea.

You may have several of predefined templates for Clinical diagnoses and for Recommendations which can be used with a click of a button to facilitate the Reviewing process.

Using a template will fill the relevant section (diagnosis or recommendation) in the report. After choosing a template you may edit the text.

Clicking on the “New” button will open the New Setting Details screen.

Please select the Type, enter the Name, and decide if Active. Then add text to the Template.

New Setting Details Screen

Existing Template Screen


4.2.2 NPIs

You can manage a list of NPI to be used later in the report by selecting Type = NPI in the search fields. Click on “Search” to see the list of all NPIs currently in the system. You can also enter text in the NPI field to limit the search results. Clicking on an existing NPI will display the NPI’s details screen, where you can edit the NPI number and save the changes.

Clicking on the “New” button and selecting Type = NPI can be used to add a new NPI to the list.

4.2.3 Signature

Interpreting Physicians can add their signature sample as a graphic file. Thus, the Interpreting Physician’s signature will be added to the report automatically. Click on the caret next to the User Name, select Set Signature and use the upload area in the screen to upload your signature file.



Interpreter Signature Screen

4.2.4 States Eligibility

Interpreting Physicians can add the list of states in which they are board certified by selecting the States Eligibility from the User menu.

4.3 Using CloudPAT as an Interpreting Physician Assistant

The Interpreting Physician Assistant is able to see the Interpreting Physician's data as read-only, and is a role meant to help the interpreting physician staff administration.

The Interpreting Physician Assistant has access to read the following:

- Studies Table
- Study Details
- Templates and NPIs Table
- Setting Details
- Pool Managers Table
- Pool Details

The Interpreting Physician Assistant has the following regulations/definitions:

The login procedure is identical to other roles. If the Interpreting Physician Assistant user is not connected to any interpreting physician – the user will not be able to log in.

The Interpreting Physician Assistant cannot create new patients or upload a study.

The Interpreting Physician Assistant receives the same emails (direct and pool related, all as bcc) as their parent but the link will lead to a read-only record.

The Interpreting Physician Assistant will see all the studies that belong to the parent – if studies are added to the parent the Interpreting Physician Assistant will see them, and if a study is taken away from the parent, The Interpreting Physician Assistant will not be able to see it.

One Interpreting Physician user can have multiple assistants connected to it.

The Interpreting Physician Assistant user role can be linked to numerous parents (Interpreting Physicians). This can be done via the "Select Parent" dropdown menu, which supports multi selection.

5 Appendices

This chapter includes the following appendices:

- [Appendix A: License Agreement](#)
- [Appendix B: Regulatory Representative](#)

Appendix A. License Agreement

License to User from Itamar Medical Ltd can be found on the Itamar Medical website:

- <https://www.itamar-medical.com/wp-content/uploads/2020/09/License-Agreement-09.2020.pdf>

Appendix B. Regulatory Representative

Itamar Medical's authorized regulatory representative is:

	Arazy Group GmbH
---	-------------------------

The Squire 12, Am Flughafen,
60549 Frankfurt am Main, Germany