

July 30, 2019

### **Approved Cleaning Products for WatchPAT Devices**

**Dear WatchPAT Customer,**

Cleaning the WatchPAT and its accessories (Snoring and Body Position sensor and Respiratory Effort and Snoring and Body Position sensor) with unapproved cleaning products may damage the device and its cables and sensors.

Below is an approved list of cleaning products to clean WatchPAT products. Please make sure that you are using one of these products to clean the device.

<b>Commercial Name</b>	<b>Manufacturer</b>
Sani-Cloth - HB	PDI
Sani-Cloth – Plus	PDI
Sani-Cloth - Bleach	PDI
Sani-Cloth – AF3	PDI
Sani-Cloth – Prime	PDI
Sani-Cloth – Super	PDI
Universal Wipes	Clinell
CaviWipes	Metrex
Clorox Bleach Wipes	Clorox

If you have questions, contact your Itamar Sales Representative or call below numbers:

#### **Inside the US**

Customer Service: **888-748-2627**

Email: [supportusa@itamar-medical.com](mailto:supportusa@itamar-medical.com)

#### **Outside the US**

Customer Service: **+972-4-6177000**

Email: [support@itamar-medical.com](mailto:support@itamar-medical.com)

Sincerely,

**Itamar Global Marketing Team**

REF: MM2198687